



## Analysis of Amines in Process of Spandex

### Analysis Condition

**Column:** HP-5MS (30m x 0.25mm x 0.25um)

**Inj. temp. :** 280°C

**Det. temp:** 280 °C

**Oven temp. (Programming):**

80°C(3min) to 150°C(1min) at 10°C/min

**Detector:** FID

**Column Flow:** 1.4ml/min

Sample	Name	RT
DEA	Diethylamine	2.001
EDA	Ethylenediamine	2.218
PDA	Propylenediamine	2.435
DMAC	Dimethylacetamide	4.213

