

APPLICATION NOTE

Analysis of Ethylene Oxide & 2-Chloroethanol

1. Analysis Condition

Acme6000 Series

Oven : 35°C(2min) –> Increasing Temp. 10°C/min –> 200°C(1min)


Column : HP-INNOWAX(30m*0.32mm*0.5um)

Carrier gas : He, 2ml/min (Split ratio 10:1)

Injector : Capillary 230°C

Detector : FID 250°C

2. Chromatogram



1. Ethylene Oxide
2. 2-Chloroethanol

