



Analysis of Epichlorohydrin

Analysis Condition

Column: HP-5 (30m x 0.32 mm i.d. x 0.25 μ m)

Injector: Temp 220 $^{\circ}$ C, Split ratio: 20:1

Detector: FID Temp: 250 $^{\circ}$ C,

Air: 300 mL/min, H₂: 35mL/min

Column flow: 2.0 mL/min

Oven: 40 $^{\circ}$ C (1min) - 10 $^{\circ}$ C/min - 100 $^{\circ}$ C

Standard: Dissolve in Chloroform



<Acme 6000GC>

