



Analysis of Acyclovir in Plasma

Analysis Condition

Pump: Younglin SP930D

Injector: Midas Autosampler (100 L loop)

Detector: LC305 Ultraflour Fluorescence

Data system: Autochro-2000

Column: Zorbax SB-18 (4.6 X 150 mm, 5 μ m)

Mobile phase: 0.02 M Phosphate buffer (pH 2.5) : ACN = 99:1

Flow rate: 1.2 mL/min

Detection: Ex: 270 nm, Em: 380 nm

